

**Surfa'safe Premium diffuse foam 750 ml**

Alcohol-free disinfectant detergent foam for cleaning and disinfecting surfaces and non-immersible, non-invasive medical devices. Microbial efficacy from 2 minutes, quick-drying and streak-free. Ready to use.

Microbiological properties: Bacteria: EN 1040, EN 13727, EN 13697: 2 mins. According to EN 13727: active on VRE, A. baumannii  $\beta$ LSE, E. coli OXA 48, E cloacae OXA 48, K. pneumoniae OXA 48, MRSA, Listeria monocytogenes, Salmonella enteritidis: 1 min. According to EN 13697: active on VRE, A. baumannii  $\beta$ LSE, E. coli OXA 48, E cloacae OXA 48, MRSA, Listeria monocytogenes, Salmonella enteritidis: 2 min. Yeasts: EN 1275, EN 13697: 5 min. EN 13624: 2 min. Mould: EN 1275, EN 13624, EN 13697: 20 min. According to EN 13624, EN 13697: Active on Aspergillus fumigatus and Tricophyton Mentagrophytes: 15 min. Viruses: EN 14476: Adenovirus, Murine Norovirus: 30 mins. According to EN 14476: active on HIV-1, PRV (HBV model virus), BVDV (HCV model virus), Vaccinia virus, HSV, Rotavirus: 1 min. According to EN 14476: active on Polyomavirus SV 40, Coronavirus: 5mn.

Dangerous - follow safety precautions. Use biocides with caution. Class IIa medical device. Read the label and product information before use.

750 ml bottle with diffuse foam dispenser.

**Date d'impression : 23/11/2024**