

Anios'Clean Excel D



- **Disinfects instruments**

Cleaning and disinfection of medical devices: surgical and medical instruments, heat-sensitive instruments and endoscopy equipment, prior to their sterilisation. Product to be diluted: 25 ml for 5 litres.

Microbiological properties: Bacteria: EN 13727, EN 14561: 5 mn - 0.5%. EN 13727 on BMR: 15 mn - 0.5%.

Mycobacteria: EN 14348, EN 14563 (M. terrae): 60 mn - 2%. Yeasts : EN 13624, EN 14562 : 15 mn - 0.5%. Viruses: PRV (HBV model virus), BVDV (HCV model virus), Herpesvirus, Vaccinia virus, HIV-1: 5 mn - 0.5%.

Dangerous - follow carefully the precautions for use. Medical device class IIb. - Before use, read the label and product information.

Available in 1 or 5 litre drums.

Variantes :

- 1 L
- 5000 ml