

## Steranios 2%



- **High-level disinfection of medical devices, surgical and medical equipment, endoscopy equipment and heat-sensitive equipment**

High-level disinfection of medical devices, surgical and medical equipment, endoscopy equipment and heat-sensitive equipment.

Microbiological properties: Bacteria: EN 1040, EN 13727, EN 14561: 5 mins. Active on *Helicobacter pylori*: 10 mins. Mycobacteria: *Mycobacterium tuberculosis* (B.K.): 5 mins. EN 14348, EN 14563 (*M. terrae*): 10 mins. Active on *Mycobacterium massiliense*: 10 mins. Yeast / Mould: EN 1275. EN 13624, EN 14562: 10 mins. Viruses: EN 14476: Poliovirus, Adenovirus, Norovirus Murine: 10 min. Active on HIV-1, HBV, BVDV (HBV model virus) and HSV: 5 min. Spores: NF T 72-230/231 *Bacillus subtilis*, *Bacillus cereus*, *Clostridium difficile*.

Dangerous - follow safety precautions. Class IIb medical device. - Read the label and product information before use.

Ready to use.

Available in 5-litre containers.