

Pic Enema bulb



- **With soft cannula : A gentle and comfortable cleansing**
- **With rigid cannula : Ideal for a simple personal hygiene**

[Pic Enema bulb](#) with rigid or soft cannula depending on the model. Suitable for auricular, intestinal and vaginal cleansing depending on the model.

Available in several models:

- 35 ml. With soft cannula. Use : auricular and intestinal cleansing.
- 70 ml. With soft cannula. Use : intestinal cleansing.
- 125 ml. With soft cannula. Use : intestinal cleansing.
- 330 ml. With rigid cannula (13 cm). Use : vaginal cleansing.

Enema spray: PVC. Non-sterile.

Before each use, wash the device with soap and water and rinse thoroughly. After use, wash with soap and water and rinse thoroughly. Allow all components to dry before storing. Before each re-use, repeat the operations specified from the beginning.

Variantes :

- 35 ml
- 70 ml
- 125 ml
- 330 ml



- 35 ml - 802067

- | | |
|-------------------|--|
| - 70 ml - 802068 | |
| - 125 ml - 802363 | |
| - 330 ml - 802364 | |

Date d'impression : 02/08/2025