

MED 528 Care Bike



- Guaranteed use 5 hours per week
- Low overhang sillMotorised resistance adjustment

Bike for rehabilitation with an aeromagnetic braking system, with motorised adjustment of 24 resistors. Inertia weight 9 kg. All comfort adjustments of the saddle and handlebars. Pulse taken by tactile sensors. Computer with LCD display with single control button. 24 preset programmes.

Dimensions : Length 95 x width 55 x height. 137 cm.

Weight: 31,1 kg.

Maximum authorised weight 120 kg.

Date d'impression : 11/04/2025