

## Universal half-tray for wheelchairs



- Right-handed half-shelf - 833038.D

- Left-handed half-shelf - 833038.G

- **Sturdy half trays that can be adapted to all types of wheelchairs**
- **2 models depending on which side of the wheelchair it is to be used.**

These half trays, made from Lexan, a material that is far more resistant than Plexiglas, can be adapted to all types of wheelchairs.

Simply slide them onto the armrest and adjust the spacing.

Dimensions : Depth 50 x Width 9/30 cm.

2 models to choose from : Left or right

**Variantes :**

- Right-handed half-shelf
- Left-handed half-shelf