

Lofoten toilet lift



- Can be installed without any specific adaptation
- Supplied with wheels: fixed or mobile installation
- Right or left hand operation

Our Lofoten toilet lift allows the user to sit down and stand up comfortably and safely from the toilet. Simply operate the control on one of the armrests and the seat will gently rise or decline. The Lofoten can be placed around your existing toilet without any fixing or renovations. It is delivered with 4 wheels (2 of which have brakes) as standard: you can choose to install it on its 4 legs (height under frame adjustable from 38.5 to 48.5 cm) or replace the legs by the wheels to make it mobile (height under frame adjustable from 42 to 45 cm).

2 models available in the Lofoten range:

The Basic model has a permanent mains connection.

The Premium model has a rechargeable battery with charge level indicator (approx. 100 cycles on a full charge), an ergonomic backrest and a removable potty for independent use.

Structure in white galvanised steel, ABS and rubber feet.

Overall dimensions: width 60.6 x depth 52.5 x height adjustable 38.5 to 48.5 cm (with feet) or 42 to 45 cm (with wheels).

Width between armrests 48 cm.

Recess : 27.5 x 23 cm.

Backrest height 38 cm (Premium version).

Maximum weight 200 kg.

All models are supplied with 4 wheels

Variantes :

- Premium
- Basic
- Right-handed
- Left-handed



- Premium - Right-handed - 811184.D
(On restock - Available from 25/06/2025)

- Premium - Left-handed - 811184.G

- Basic - Right-handed - 811183.D

- Basic - Left-handed - 811183.G