

## Moovy Electric toilet/shower and transfer chair

- The electric lift version of the [Moovy chair](#)

The electric lift version of the [Moovy chair](#)

Like the [standard Moovy chair](#), it allows you to take charge easily, without effort for the caregiver and without risk of pain for the patient. It will allow you to transfer your patients from bed to shower or toilet.

Simply position the seat of the Moovy at the height of the patient's buttocks, by changing the height using the electric jacks. The seat can be raised and lowered directly on the Moovy (without remote control). The seat must then be opened, closed around the patient and secured with the double safety system (cleat and belt buckle).

It has 4 wheels with brakes, a soft safety belt and soft foam back and armrests.

It comes with a comfort cushion, and a removable potty for independent use.

It also has folding handles for better access to the patient.

It can be used in the bathroom shower as long as the battery is removed before use. Please refer to the instructions and adhere to the safety precautions and operating guidelines for further information.

Overall dimensions: width 62.5 x depth 78 cm (floor) or 82 cm (top).

Seat dimensions: width 47 x depth 43 x height. 44.5 to 64 cm. Back height 40.5 cm.

Width between armrests 45.5 cm.

Weight 35 kg.

Maximum weight 150 kg.



- 822013