

**Groenland folding toilet frame**



● **Folds flat for easy storage**

Our Groenland toilet frame allows you to sit at the right height and stand up easily from the toilet. Can be used on a toilet bowl or on its own with the supplied potty. It has armrests and a backrest. The entire unit folds flat. Supplied with a potty and removable lid. Equipped with a guide to position it over the toilet.


Overall dimensions: width 56.5 x depth 53.5 x height 61 to 76.3 cm.

Weight : 5.2 kg.

Seat dimensions: width 33.5 x depth 37 x height. 39.5 to 55 cm.

Width between armrests : 46 cm.

Maximum user weight : 135kg.



Comment sizer le cadre de toilette ? 1/2



Comment plier le cadre de toilette ? 2/2



Comment s'installer sur toilette ? 3/2



- 811187