


Groenland folding toilet frame



• **Flat folding for easy storage**

The Groenland toilet frame allows a user to sit at the right height to get up easily from the toilet. It can be used over a toilet bowl or on its own with the supplied potty. It is equipped with armrests for support and a backrest. The entire unit is flat folding. Supplied with a removable potty with a lid with a removable from the rear of the chair. Equipped with a guide to position it over the toilet.

Overall dimensions : width 55 x depth 51 x height 59 to 78 cm.
 Seat dimensions : width 33.5 x depth 37 x height 39.5 to 55 cm.
 Width between armrests : 46 cm.
 Weight : 5.2 kg.
 Maximum user weight : 135 kg.



Comment s'ajuste le cadre de toilette ? 1/2

Comment placer le cadre de toilette ? 2/2

Comment s'ajuste le cadre de toilette ? 3/2



- 811187