

## Air comfort mattress overlay



This mattress overlay provides perpetual movement on the bed. This movement allows people who are permanently bedridden to change position. The mattress is made up of two separate rows of cells spread across its entire surface.

The movements are divided into three phases :

Phase 1 : line 1 inflates to maximum and line 2 deflates.

Phase 2 : both lines inflate to maximum.

Phase 3 : line 1 deflates and line 2 inflates to maximum.

Variable pressure compressor.

Pressure adjustable from 50 to 110 mmHg.

Number of litres/minute : 5 to 6.

Duration of a cycle (3 phases) : 12 minutes.

Mattress dimensions : Length 200 x width 90 x depth 7 cm.



- 823024

**Date d'impression : 10/04/2026**