

Retractable bed rail



- Helps to prevent restless people from falling during sleep
- Operable by the user

This bed rail can be operated by the user. It is fixed to the bed base and helps to prevent restless people from falling out of bed while they are asleep. It can be folded away with 2 easy-to-grip clips and is available in 3 different lengths with just one clip (66, 86 and 107 cm). Adaptable on a box spring or slatted bed base. Delivered with a storage net. The retractable bed barrier works with most types of traditional beds with mattresses of thickness between 30 and 40 cm.

Total height: 58.4 cm.

Height of the retractable part: 48 cm. Height of the barrier: 33 cm.

Weight 6.35 kg.

Maximum supported weight: 90 kg

Date d'impression: 21/11/2024