

Capiluve shampoo basin for bedridden patients



- **Shampoo basin for easy and effortless shampooing of bedridden patients**

The Capiluve makes it very easy to wash the hair of bedridden patients. Simply slide it under the patient's shoulders so that it stays firmly in place, and position the neck in the basin. The wide basin allows for freedom of movement. Water does not stagnate, but flows into the Capiluve's internal reservoir, which can hold up to 8 litres. When the shampoo is finished, the entire unit can easily be transported without any water leakage and emptied using the drain plug. Effortless shampooing made easy.

Dimensions : 60 x 60 x 16 cm.

Weight : 2.7 kg.

Material: Plastic polyethylene



- 823059

Date d'impression : 06/02/2026