

## Heated lumbar belt



- Self-heating - 835180
- Hot Cold - 835179
- Thermo-cryo - 835181

- **3 models available**

Very easy to use, simply attach the lumbar belt to the waist to benefit from the heat on the lumbar region. This belt can also be used wrapped around a thigh or an arm. Length: 120 cm.

3 models available:

- Hot cold model with microwave heating
- Auto model with automatic heating
- Thermo-cryo model with bowl heating.

**Variantes :**

- Self-heating
- Hot Cold
- Thermo-cryo