

Home press-therapy Device



- Standard - 838005

- Pro - 838013

Revitalising massage by pressure and lymphatic drainage. Improves venous circulation by restoring venous function. Promotes blood circulation and relieves heavy and tired legs. Prevents varicose veins and spider veins. Massage by air pressure and inflation and deflation of the cushions. Intensity of the massage can be adjusted continuously and alternately on the leg pads. Timer included: 10/20/30 min.

Storage bag included.

Battery operated (3 AAA batteries included) or mains operated.

Legrest dimensions: 61 x 67 cm.

Warranty: 3 years.

2 models available:

- **Standard : leggings only.**
- **Pro : full legs with carry bag (Pro model)**

**Variantes :**

- Standard
- Pro