

Sharpsafe collector



- **11 sizes to choose from**

Portable containers ranging from 0.2L to 0.6L for emergency services.

0.2 L container (needle disconnection + injection pen)
 0.3 L and 0.45 L (needle disconnection)
 0.6 L (full opening without disconnection).

Containers from 1 L to 13 L for sharp, cutting and piercing objects. 1 L and 13 L containers (without flap), from 2 L to 9 L (lid with flap). Wide opening 55x35 (1 L), 55x50 (2, 3, 5 and 7 L), 210x65 (9 and 13 L).

0.2 L : Dimensions : 146 x 79 x 37 mm.

0.3 L : Dimensions : 130 x 105 x 50 mm.

0.45 L : Dimensions : 170 x 105 x 50 mm.

0.6 L : Dimensions : 224 x 108 x 52 mm.

1 L : Dimensions : 187 x 90 x 108 mm.

2 L : Dimensions : 177 x 197 x 126 mm.

3 L : Dimensions : 203 x 197 x 126 mm.

5 L : Dimensions : 260 x 170 x 205 mm.

7 L : Dimensions : 380 x 175 x 175 mm.

9 L : Dimensions : 291 x 292 x 195 mm.

13 L : Dimensions : 353 x 292 x 195 mm.

Variantes :

- 1 L
- 13 L
- 0,2 l
- 2 L
- 0,3 l
- 3 L
- 4 L
- 0,45 l
- 5 L
- 0,6 l
- 7 L
- 9 L