

## Nosy Be toilet frame



- **Folding armrests**
- **Facilitates transfers**

The Nosy Be toilet frame can be used independently with the provided potty, or over an existing toilet. It has padded armrests that can be folded away for easy lateral transfers (folding by means of a simple push button).

Overall dimensions : width 55 x depth 52 x height 72 to 86.5 cm.


Seat dimensions : width 35.5 x depth 43 x height adjustable from 45 to 57 cm in 7 positions.

Width between armrests : 46.5 cm.

Backrest height : 25 cm.

Weight : 5kg

Maximum user weight : 136 kg



- 811011

Date d'impression : 06/02/2026