

## **Nosy Be toilet frame**



armrests that can be folded away for easy lateral transfers (folding by means of a simple push button). Overall dimensions: width 55 x depth 52 x height 72 to

The Nosy Be toilet frame can be used independently with the provided potty, or over an existing toilet. It has padded

Seat dimensions: width 35.5 x depth 43 x height adjustable

from 45 to 57 cm in 7 positions. Width between armrests: 46.5 cm. Backrest height: 25 cm.

 Folding armrests • Facilitates transfers

Weight: 5kg

86.5 cm.

Maximum user weight: 136 kg













**Date d'impression :** 02/08/2025