

Majorque home patient lift



- 825012

- **Supplied with a strap size M**

Manual foot spreader. 4 castors: 2 x 75 mm at the front and 2 x 100 mm with brakes at the back. Equipped with a comfortable double padded hand rest. Remote control with stand. Steel frame. DC 24V 5 Ah battery. DC 24V 8000N motor. Full charge 5 hours. Each patient lift is delivered with a size M strap as standard.

Dimensions:

Useful length of legs 68 cm, total length 95 cm.

Width of legs in straight position 54 cm, in spread position 88 cm.

Height of legs 95 mm.

Height of the bottom of the flail minimum 69 cm, maximum 171.5 cm.

Maximum operating weight: 150 kg

Flail dimensions :

Length 49.2 cm x width 41 cm.

Width between strap supports: 39 cm.

Weight 40 kg.