

HEBIX LIFTING AID CHAIR



- Dismountable chair - 821040

- Transport bag - 821106

- Chest strap - 821109

- Folding chair - 821041

(On restock - Available from 18/04/2025)

- **A simple and innovative solution to gently lift people off the floor**
- **Ideal to relieve the carers**

In a few simple steps, the Hebix can be opened and placed on the person lying on the floor. The decisive factor here is that the force required is only about a quarter of the body weight of the person to be lifted: for a person weighing 80 kg, it is only about 20 kg. By simply lifting the Hebix using the handles, you can gently and painlessly move the person lying down into a sitting position. The person can remain in this position, stand or be moved, nothing could be easier !

2 models available:

- A 3-part dismantled model. Weight 7 kg.

- Foldable model, more suitable for emergency services. Weight 10.5 kg.

Dimensions : Length 75 x width 45 cm.

Resistance : 100 kg.

The optional carry bag is suitable for the folding model only. You can opt for a chest strap to secure the patient during lifting.

Resistance : 100 kg

Variantes :

- Dismountable chair
- Transport bag
- Chest strap
- Folding chair