

## **Pressotherapy I-Press II**





- M Device 838155 (On restock - Available from 27/11/2024)
- M Boot expander 838157.M
- L Device 838156
- L Boot expander 838157.L
- One size Right sleeve 838158.D
- One size Left sleeve 838158.G
- One size Tubing 838177

The I-Press II pressure therapy device is ideal for relieving circulatory system disorders to improve peripheral blood circulation. Thanks to its pumping action, it promotes venous return.

### **Technical specifications:**

Max nominal pressure 200mmHg +-20%, max absorption 100 mA. Therapy time adjustable from 0 to 30 minutes. Mains powered. 2 year warranty. The device is delivered with a pair of boots and 2 massage insoles.

#### 2 different sizes:

- Size M:  $\max$  thigh circumference 60 cm,  $\max$  ankle circumference 33 cm.
- Size L: max thigh circumference 72 cm, max ankle circumference 46 cm.

# The options of the I-Press II pressure therapy device

- The boot expanders allow the use of the pressure therapy device when the existing boots are not wide enough. Gain 13 cm of thigh circumference and 5 cm of knee circumference. Available in sizes M and L
- The sleeves allow to work the arms (arm circumference 46 or 59 cm, armpit to wrist length 87 cm). Available for the left and right arm. One size fits all.
- The belt allows you to work the abdominal and lumbar muscles (waist circumference 117 to 132 cm). One size fits all.

With the purchase of the sleeves and the belt, you will necessarily need a tubing to connect them to the I-Press II pressotherapy device.

#### Variantes:

- M
- L
- One size
- Device
- Belt
- Boot expander
- Right sleeve
- Left sleeve
- Tubing

Date d'impression: 21/11/2024