

Retractable bed rail with stand



- 823157

(On restock - Available from 09/04/2025)

- **Very safe due to its 2 legs on the ground**
- **Comes with a safety strap and buckle**

This practical and reassuring barrier can be lowered simply by lifting the foam handle upwards. The length of the barrier can be adjusted in 3 lengths for more or less protection: 54.5 - 64.5 - 74.5 cm. Partially upholstered top support for easy removal from bed.

Height of the upper support in relation to the underside of the mattress 55 cm.

Dimensions of the part that fits between the mattress and the base: 57.5 x 55.5 cm.

Height of legs from floor to underside of mattress: min. 22.2 /max. 34.9 cm, continuously adjustable.

Strap with buckle, length 384 cm.

Weight : 4.3 kg.

Date d'impression : 23/02/2025