

## **Retractable bed barrier with stand**



- Very secure due to its two sturdy legs
- Comes with a safety strap

This practical and secure barrier can easily be removed by lifting the foam handle upwards. The length of the barrier can be adjusted to three different lengths to provide varying degrees of protection: 54.5 cm, 64.5 cm and 74.5 cm. The top rail is partly padded to make it easier to get out of bed.

Height of the upper support from the bottom of the mattress :  $55\ \text{cm}$ 

Dimensions of the part that fits between the mattress and the bed base :  $57.5 \times 55.5 \text{ cm}$ .

Height of the legs from the floor to the bottom of the mattress: min. 22.2 / max. 34.9 cm, continuously adjustable.

Strap with buckle length: 384 cm.

Weight: 4.3 kg.



 $\textbf{Date d'impression:}\ 15/12/2025$