

## **GelPro hot-cold compresses shoulder**



- Revolutionary stretch material!
- Enables efficient cold or hot treatment

The new GelPro range offers a complete thermotherapy treatment. By placing the compress in a freezer, you can relieve minor trauma and oedema. By placing the compress in the microwave, you can soften and relieve muscular or body pains.

The compress fits closely over the area to be treated because the material used is extensible, allowing it to be adapted to most body shapes.

Blue colour.

Half-belt, half-tube model that slips over the arm to cover the entire shoulder. Stretch velcro strap secures under opposite arm.

Dimensions without extension: 29 x 35 cm. Strap length without extension 76 cm.

Compresses available in 5 other models.

- 835192 (On restock - Available from 19/03/2025)

 $\textbf{Date d'impression:}\ 23/02/2025$