

## Madagascar transfer chair



- 822275

- **The new generation of transfer devices**
- **Allows the user to shower without further transfer**

The Madagascar Transfer Chair allows for quick, effortless, and safe transfer of a person without causing discomfort. From a seated position, simply slide the soft seat under the user and secure the four side straps. After positioning the safety backrest, lift the user electrically. The battery-powered movements are smooth and steady. Adjustable leg spacing provides additional safety during use. The Madagascar allows the user to be transferred from one seat to another, and also enables hygiene care without additional transfers. Quick release is possible by pressing the four side straps. Equipped with 2 wheels with brakes. Seat height adjustable via the adjustable seat straps.

Overall dimensions: W 60.5 × D 87 cm.

Leg spacing: 60.5–91 cm.

Usable width within the frame: 51 cm.

Handle height: 86.5–122 cm.

Removable backrest: W 33 × H 20 cm.

Weight: 33.5 kg.

Operates on a rechargeable 2Ah 44.4Wh battery with a DC 24V motor.

Maximum user weight: 150 kg.